

# Rapicide® OPA/28

HIGH-LEVEL DISINFECTANT

Works faster.  
Lasts longer.





# When safety counts and minutes matter, choose RAPICIDE® OPA/28

Ensuring a safe and healthy environment for patients and staff just got easier. RAPICIDE OPA/28 provides a superior way to reprocess heat-sensitive, semi-critical medical devices that are unsuitable for sterilization. It features the fastest disinfection time, twice the reuse period of other OPA brands and guaranteed materials compatibility, which allows for the ultimate combination of safety, convenience and value.



## The RAPICIDE OPA/28 Advantage® **BY THE NUMBERS**

**5** MINUTE AER disinfection time at 25°C

**10** MINUTE manual disinfection time at 20°C minimum

**84** reuse cycles<sup>1</sup>


**28** DAY reuse period

**2** final rinses needed in MEDIVATORS® AER

**3** final rinses needed for manual applications

## **WORKS FASTER**

Take quicker action against microorganisms



With RAPICIDE® OPA/28, you have a powerful, ready-to-use disinfectant carrying an active ingredient of 0.575% ortho-phthalaldehyde (OPA). High-level disinfect devices manually in just 10 minutes at room temperature (20°C minimum), to achieve bactericidal, fungicidal and tuberculocidal kills.

Leverage proprietary surfactant technology for superior results

RAPICIDE OPA/28's proprietary surfactant technology aids disinfectant penetration through cell walls of microorganisms — improving disinfection efficacy relative to other OPA brands. Studies<sup>2</sup> indicate it effectively inactivates:

- MRSA
- VRE
- TB
- HIV
- Hepatitis Viruses



Surfactant technology at work<sup>3</sup>

# LASTS LONGER

## Use for up to 28 days

It's important for healthcare providers to manage budgets carefully while still delivering high-quality patient care. RAPICIDE® OPA/28 High-Level Disinfectant has been tested to provide 84 reuse cycles<sup>1</sup> along with a reuse period of up to 28 days — twice the reuse period of other OPA brands. This allows the use of one bottle versus two<sup>4</sup> over a 28-day period, creating potential cost-saving opportunities and protecting your workers by limiting exposure to disinfectants.




No mixing or activation is needed.



# DOES MORE

## Get more disinfection cycles in the same amount of time

RAPICIDE OPA/28's faster manual disinfection time provides one extra disinfection cycle for every five instruments reprocessed, allowing you to increase your operational efficiencies and reduce downtime of medical instruments.



## To order

Contact your distributor or call one of these Cantel Medical companies:


SPSmedical: (800) 722-1529

Crosstex: (888) 276-7783

## Test Strips Provide MRC Verification

RAPICIDE® OPA/28 Test Strips are designed for specific use with our disinfectant, and make measuring the concentration:

- **User friendly**—strips feature a large, easy-to-read indicator pad
- **Intuitive**—color chart on bottle provides a pass/fail reference
- **Fast**—3-second dip time and 90-second test strip read time



<sup>1</sup>MEDIVATORS Internal Testing Per EPA Reuse Test Protocol Specification

<sup>2</sup>Test data pertaining to all microorganisms is available upon request

<sup>3</sup>AER use may generate minimal suds, which do not compromise worker or patient safety. Follow proper disinfectant chemical replacement procedure and equipment maintenance for best disinfection and rinse results.

<sup>4</sup>Provided that MRC meets requirements and product is used for full 28-day period

# Trust in proven quality from the MEDIVATORS® RAPICIDE® disinfectant family

You can trust that RAPICIDE OPA/28 matches your high standards for patient care and staff safety. It's manufactured by MEDIVATORS, a proven market leader with the highest customer satisfaction rating in the endoscope reprocessing industry.<sup>1</sup> You can use our disinfectant manually or seamlessly with selected MEDIVATORS reprocessors or other compatible reprocessors.



## RAPICIDE® OPA/28 High-Level Disinfectant

<b>Part Number</b>	<b>Description</b>
ML02-0127	RAPICIDE OPA/28 High-Level Disinfectant (4 x 1-gallon bottles/case)
ML02-0137	RAPICIDE OPA/28 Test Strips (2 x bottles/package; 50 strips/bottle)

### To order

Contact your preferred distributor or a representative from one of the Cantel Medical companies below.

#### Distributed by:



Crosstex International, Inc. 10 Ranick Road Hauppauge, NY 11788 (888) 276-7783 www.crosstex.com	SPSmedical, a Crosstex Company 6789 W. Henrietta Road Rush, NY 14543 Toll Free: (800) 722-1529 www.spsmedical.com
---	---

#### Manufactured in the USA by:



MEDIVATORS Inc.  
14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Toll Free: (800) 444-4729  
www.medivators.com

<sup>1</sup>MD Buyline Intelligence Report, December 2012

MEDIVATORS® and RAPICIDE® are registered trademarks of Medivators Inc.  
© 2015 Medivators Inc. P/N 50098-075C

