

Black Line Rubber Dam Clamps

An Exceptional Restoration Begins with Excellent Isolation



Instrumentation and
Procedural Products

Our **Rubber Dam Clamps**,
made in the USA from proprietary stainless steel,
provide full-mouth coverage.

Anterior Rubber Dam Clamps

<p>B4X RDCMB4X</p> <p>Anterior incisors and canines.</p>	<p>106RX RDCM106RX</p> <p>Anterior region.</p>	<p>106LX RDCM106LX</p> <p>Anterior region.</p>
<p>6X RDCM6X</p> <p>Slightly offset jaws for labial caries on anterior teeth.</p>	<p>9X RDCM9X</p> <p>Offset flat jaws and stiff bow for anterior teeth.</p>	<p>9SX RDCM9SX</p> <p>Offset flat serrated jaws for anterior teeth.</p>
<p>00X RDCM00X</p> <p>Flat jaws and high bow for small lower canines and first premolars.</p>	<p>212X RDCM212X</p> <p>Flat jaws for labial caries on anterior teeth.</p>	<p>212SAX RDCM212SAX</p> <p>Slightly offset staggered jaws for odd shaped anteriors.</p>



Satin Steel™ Alloy
which **resists
corrosion while
aiding in flexibility
and memory**

A dark coating allowing
for **improved intraoral
contrast and visual acuity**

Matte finish for
**decreased light
reflection and glare**

Molar Rubber Dam Clamps (General)

<p>8AX RDCM8AX</p> <p>Deeply festooned four-point jaws for molars.</p>	<p>8ADX RDCM8ADX</p> <p>Deeply festooned jaws and extended bow for distal access to molars.</p>	<p>W8AX RDCMW8AX</p> <p>Slightly festooned jaws for molars.</p>	<p>10X RDCM10X</p> <p>Flat jaws designed for small upper left and lower right molars; paired with RDCM11.</p>	<p>11X RDCM11X</p> <p>Flat jaws designed for small upper right and lower left molars; paired with RDCM10.</p>	<p>14X RDCM14X</p> <p>Deeply festooned jaws for partially erupted or irregularly shaped molars.</p>	<p>14AX RDCM14AX</p> <p>Deeply festooned jaws for partially erupted or irregularly shaped molars; larger than RDCM14.</p>
<p>W14AX RDCMW14AX</p> <p>Deeply festooned jaws for partially erupted or irregularly shaped molars.</p>	<p>18X RDCM18X</p> <p>Flat jaws and stiff bow for molars.</p>	<p>24X RDCM24X</p> <p>Flat jaws for upper right and lower left molars; paired with RDCM25.</p>	<p>25X RDCM25X</p> <p>Flat jaws for upper left and lower right molars; paired with RDCM24.</p>	<p>27X RDCM27X</p> <p>Flat jaws for large molars.</p>	<p>56X RDCM56X</p> <p>Slightly festooned jaws for large molars.</p>	<p>56SX RDCM56SX</p> <p>Slightly festooned serrated jaws for large partially erupted or irregularly shaped molars.</p>

Black Line Rubber Dam Clamps



Instrumentation and
Procedural Products

Complementary Products



Scan or click to
Request Your
Dental Dam
Sample

We also offer powder-free Dental Dams,
available in latex and non-latex options.



Rubber Dam
Forceps
RDF

Rubber Dam
Punch-Ainsworth
RDPA


Rubber Dam
Punch
RDP




Rubber Dam
Frame, Adult
RDAF6

Rubber Dam
Frame, Child
RDCF5


Premolar Rubber Dam Clamps



Upper Molar Rubber Dam Clamps



Lower Molar Rubber Dam Clamps



To learn more about our Black Line Rubber Dam Clamps, visit: HuFriedyGroup.com/BlackLineRDC

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com
 All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.
 Marks not registered in all jurisdictions.
 ©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HF-5696 /1023

HuFriedyGroup
The Best In Practice