

Power Plus Ultrasonic Inserts

Reach Further for Pocket Debridement



New **Power Plus Inserts** feature **Power** for effective removal of tenacious deposits and heavy extrinsic stain **PLUS** extended length to reach deep into periodontal pockets. With a 25% longer tip design,* a variety of tip sizes, and a focus on clinician comfort, Power Plus Inserts offer an advanced solution for debridement of tenacious deposits.



Power Plus Standard Beveled Swivel Insert
UI30PPSB



Power Plus Thin Beveled Swivel Insert
UI30PPTB



Power Plus Standard Canonical Swivel Insert
UI30PPSC











Power Plus Thin Canonical Swivel Insert
UI30PPTC

FEATURES

- **2+ mm longer tip*** for extended reach into periodontal pockets
- **Tip options** to accommodate diverse pocket anatomy
 - Thinner tip for tenacious deposits in narrow pockets
 - Standard size tip to reach heavy debris in deeper pockets
- **Tip shape** similar to periodontal probe for reach, familiarity and ease of use
- **Special blend of stainless steel** for durability and tip longevity
- **Improved clinician comfort**
 - Popular Swivel™ handle for excellent control and maneuverability
 - Ergonomic wide diameter handle for reduced finger pinching
 - Textured silicone grip for secure grasp
 - Uniform diameter handle for ease of movement when shifting grip location

Power Plus Ultrasonic Inserts



INSERT TIP	TIP DESCRIPTION	RECOMMENDATION	USAGE TIPS	POWER SETTING	HANDLE TYPES	WATER DELIVERY	HANDLE COLOR	GRIP	PART CODES
Power Plus Standard Size Conical Tip 	Long, standard diameter shank extending to a tapered tip	<ul style="list-style-type: none"> Moderate to heavy deposits; supragingival; wide, deep pockets 2 mm longer than #10 insert 	Deposit removal in deep pockets; areas where gingival tissues is distended For use on moderate to heavy deposits	Low to High	Swivel	Base Flow	Ocean Blue with black cone	Silicone	UI30PPSC
	Tip cross section: 								
Power Plus Thin Conical Tip 	Long, thin diameter shank extending to a tapered tip	<ul style="list-style-type: none"> Light to moderate deposits; narrow, deep vertical pockets; areas mesial to retromolar pad; crowded areas; general supragingival use 2 mm longer than #10 insert 	Excellent choice for a variety of clinical situations. Thin shank allows effective access to narrow, deep vertical pockets as well as developmental concavities and tight areas between misaligned teeth. Ideal for localized periodontal defects and periodontal maintenance therapy.	Low to Medium	Swivel	Base Flow	Ocean Blue with black cone	Silicone	UI30PPTC
	Tip cross section: 								
Power Plus Standard Size Beveled Tip 	Long, standard diameter shank extending to a tapered, beveled tip	<ul style="list-style-type: none"> Moderate to heavy deposits; tenacious deposits; supragingival; wide, deep pockets 2 mm longer than #10 insert 	Ideal for removal of tenacious calculus in deep subgingival pockets particularly where gingival tissue is distended. Beveled edges allow for concentrated power to fracture tough deposits.	Low to High	Swivel	Base Flow	Ocean Blue with black cone	Silicone	UI30PPSB
	Tip cross section: 								
Power Plus Thin Beveled Tip 	Long, thin diameter shank extending to a tapered, beveled tip	<ul style="list-style-type: none"> Light to moderate deposits; narrow, deep vertical pockets; crowded areas 2 mm longer than #10 insert 	Great for narrow, deep vertical pockets with light to moderate tenacious calculus. Beveled edges concentrate energy on corners for enhanced deposit removal.	Low to Medium	Swivel	Base Flow	Ocean Blue with a cone	Silicone	UI30PPTB
	Tip cross section: 								

For more information visit [HuFriedyGroup.com/Inserts](https://www.HuFriedyGroup.com/Inserts)

*Compared to Hu-Friedy #10 Universal

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | [HuFriedyGroup.com](https://www.HuFriedyGroup.com)

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.

©2022 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-308/0922