

Introducing the Hu-Friedy

## BLACK LINE SCISSORS COLLECTION



Designed for performance, the Hu-Friedy Black Line is engineered to optimize clinical outcomes by delivering efficiency throughout the entire perio and surgical procedure. From reducing light reflection to enhancing visibility at the surgical site, the Black Line was designed with the clinician in mind. In addition to the surgical hand instrument portfolio, the Black Line now includes 7 Super-Cut scissor designs delivering superior performance through enhanced precision and innovation.

### THE BLACK LINE SCISSORS FEATURES:

- Super-Cut design featuring a **razor sharp blade and serrated edge for reduced tissue slippage**
- Performance engineered coating for **optimal edge retention**
- Distinct black finish for **enriched contrast and visual acuity at the surgical site and underlying tissue**
- **Reduced light reflection** afforded by a matte black finish
- Gold screw for easy identification of **Super-Cut design**

GOLD SCREW FOR  
EASY SUPER-CUT  
IDENTIFICATION

SUPER-CUT DESIGN  
FEATURING A RAZOR SHARP  
BLADE AND SERRATED  
EDGE FOR REDUCED TISSUE  
SLIPPAGE

PERFORMANCE  
ENGINEERED COATING  
FOR REDUCED LIGHT  
REFLECTION AND  
ENHANCED CONTRAST

How the best perform

**Hu-Friedy**

# BLACK LINE SCISSORS



LaGrange  
Double-Curved,  
Super-Cut

| S14SCX



Metzenbaum  
Super-Cut

| S19SCX



Iris,  
Straight/Delicate,  
Super-Cut

| S17SCX



Iris,  
Curved/Delicate,  
Super-Cut

| S18SCX



Goldman-Fox,  
Straight,  
Super-Cut

| S20SCX




Goldman-Fox,  
Curved,  
Super-Cut

| S16SCX



Castro Microsurgical,  
Super-Cut

| SPVX



Visit us online at [Hu-Friedy.com](http://Hu-Friedy.com)

©2014 Hu-Friedy Mfg. Co., LLC. All rights reserved. Images are to scale. M9090/0814

**Hu-Friedy**