

HuFriedyGroup Key Opinion Leader **Dr. Sascha A. Jovanovic**



The Implant Surgery & Tissue Regeneration Kit



Focus on:
**Periodontal & Implant
Surgery**



Country: **US - California**



Favourite Instrument: **PJO/16
JO 1 - Jovanovic 4 mm Straight/
Backaction Chisel**



Website: **www.gidedental.com**



Inspired by:
Prof. Hubertus Spiekermann



University:
**UCLA, LOMA LINDA
& ACTA**



VISIT THE
KOL COMMUNITY



SJOVANKIT
Full Jovanovic Implant
Surgery & Tissue
Regeneration Kit

Visit us online at **HuFriedyGroup.com**




HuFriedyGroup instruments are the best high quality, precise hand instruments in the industry, and a natural extension of my hands.

Dr. Sascha A. Jovanovic



THE IMPLANT SURGERY & TISSUE REGENERATION KIT

The Jovanovic Implant & Regeneration Surgical Kit includes exceptionally crafted instruments specifically designed for advanced & minimal invasive periodontal and implant surgery. This sophisticated kit includes precision-engineered instruments fashioned from premium materials. Ergonomics plays a critical role in a clinician's long-term success which led to the incorporation of Black Line instruments. The dark coating minimizes glare, is visible when used subgingivally, and adheres less to surgical materials. The use of IMS™ Cassettes ensures easy access to instruments, keeping everything needed for procedures within reach.



THE IMPLANT SURGERY & TISSUE REGENERATION KIT



18

10-130-05E

Scalpel Handle
#5E European Style
Round, Matte Finish



19

SHDPV

Scalpel Handle Steel
Microsurgical



12

NHDPV

Diamond Dusted
Castroviejo
Microsurgical Needle
Holder, Straight, 18 cm



13

TPSLCOSM

Micro Tissue Plier,
Straight, 17.78 cm



11

PFIWDS1MKX

Microsurgical Periosteal
Black Line



17

MIR5HDB (2x)

#5 HD **Black Line**
Single-Sided Mirror

16

MH6B (2x)

#6 **Black Line** Cone
Socket Mirror Handle



9

CRMX

University of Minnesota
Retractor, **Black Line**

BLACK LINE

Black Line is designed to optimize clinical outcomes by delivering efficiency throughout the entire perio and surgical procedure. Black Line features a performance engineered coating for optimal edge retention, reduced light reflection afforded by the black matte finish and enhanced contrast and visual acuity at the surgical site and underlying tissue.



22

PESIM16

Sinus Lift Instrument
SE #6



4

CL86X

#86 Lucas Spoon
Shape Surgical
Curette, DE
Black Line



3

SKRA36

#KRA-3 Kramer
Periodontal Surgical
Curette DE,
#6 Handle



23

SGCSUB-06

Sub-0 Vision
Curvette #6



24

SC4R/4L9E2

Columbia 4R/4L
Curette #9



10

CRMS

Simion Cheek
Retractor



14

S14SCX

14 LaGrange Scissors,
Curved, Super-Cut
11.5 cm (4.5"), 1 serrated blade
Black Line



28

S20SCX


Goldman-Fox Scissors, Straight, Super-Cut
13.5 cm (5.25"), 1 serrated blade
Black Line

SUPER-CUT

Cutting blade honed/
serrated to "razor"
sharpness with opposing
blade finished with fine
serrations to avoid tissue
slippage.

THE IMPLANT SURGERY & TISSUE REGENERATION KIT

INSTRUMENTS DESIGNED IN COLLABORATION WITH DR. SASCHA A. JOVANOVIC



PJO/16

JO 1 - Jovanovic
4 mm Straight/
Backaction Chisel
#6 handle

PJO/2M6

JO 2 - Mini Jovanovic
2 mm Straight/
Backaction Chisel
#6 handle

PJO/3M6


JO 3 - Mini Jovanovic
3 mm Straight/
Backaction Chisel
#6 handle

PLGJO/4X

Jovanovic Bone Graft
Condenser DE
3 mm, Black Line

PLGJO/56

Jovanovic Bone Graft
Carrier/Condenser
#6 handle



IMPJO66

Jovanovic Sinus
Lift Instrument
#6 Handle

IMPJO76

Jovanovic Sinus
Lift Instrument #7,
#6 handle



TP5090

DeBakey Perma Sharp™
Tissue Plier with teeth,
15.24 cm



TP5070

General Perma Sharp™
Tissue Plier no teeth,
15.24 cm



NH5024R

Perma Sharp™ Castroviejo
Needle Holder, Round, 18 cm



PCPUNC15CH

Quilix SE UNC Probe, 1-15 Markings, #Smooth Handle



PPR3X

PR-3 Prichard
Periosteal, **Black Line**

PERMA SHARP™ INSTRUMENTS
Featuring tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.

CREATE YOUR CUSTOM KIT AND CHOOSE THE CASSETTE WITH YOUR FAVORITE COLOR.


Just add the color coding number to the last part of the partcode (replace " _ " with the number)

| | |
|------------|----|
| Red | 0 |
| Grey | 1 |
| White | 2 |
| Orange | 3 |
| Purple | 4 |
| Yellow | 5 |
| Lavender | 6 |
| Black | 7 |
| Blue | 8 |
| Green | 9 |
| Ocean Blue | 11 |

WHY YOU SHOULD USE AN INSTRUMENT MANAGEMENT SYSTEM (IMS):

By standardising your instruments with IMS System, additional time and revenue is generated while extending the useful life of your HuFriedyGroup instruments. Dental professionals state that by implementing IMS™ Cassettes, they save 5 to 10 minutes per procedure allowing the staff to spend more time on value-added activities*. In addition, safer transportation of soiled instruments reduces risk of injury and exposure to blood borne pathogens.

*When compared to single instrument decontamination. Based on market survey results.



29 IMN4203

IMS™ Infinity Series™ **20 Instrument Cassette, Orange**
Dimensions: 203 mm x 280 mm x 32 mm

SUGGESTED INSTRUMENTS (not included in the kit)



CM2/4X
2/4 Molt DE Curette **Black Line**



SG11/12R9E2
Gracey 11/12 Curette
Rigid #9



SG13/14R9E2
Gracey 13/14 Curette
Rigid #9



SG5/6R9E2
Gracey 5/6 Curette
Rigid #9




SG7/8R9E2
Gracey 7/8 Curette
Rigid #9

SUGGESTED INSTRUMENTS (not included in the kit)



RBOE
Boehler Rongeurs
double hinged, 16 cm



H4
#4 Mosquito Hemostat,
Straight 12 cm



DID YOU KNOW?

All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilized in the open position.**

SUGGESTED BLADES AND SUTURES (not included in the Kit)



PSNPTFE1
3-0 PTFE, Undyed Reverse Cutting 19 mm
3/8 circle Reverse Cutting 45 cm



40-815
15 Sterile Scalpel Blades, QTY: 100/box,
Surgical Grade Stainless Steel



40-815C
15C Sterile Scalpel Blades, QTY: 100/box,
Surgical Grade Stainless Steel



MB69
Full Radius Round Scalpel Blades, Sterile,
QTY: 12/box, Surgical Grade Stainless Steel



MB67
Curved Edge Mini Scalpel Blades, Sterile,
QTY: 12/box, Surgical Grade Stainless Steel

Knowledge is power



Subscribe to the **gIDE** Video Streaming Platform



\$399.00
per Year

\$49.00
per Month

Evidence Based **Video Education** at its best.



Subscribe here

- 550+ Clinical Videos on Web and App
- New Dental Videos Added Weekly
- Implants, Periodontics, Esthetics & Digital
- For General Dentists & Specialists
- Reduced Pricing on Classroom Training

The gIDE institute is the premiere provider of online dental education with members and online users in over 100 countries and a world renowned faculty of clinicians and researchers from US, Europe and Australasia. gIDE has mastered the technique for delivering high quality digital media - Lecture-On-Demand, Clinical Videos, Expert interviews, Live Broadcasts over the internet to dental professionals, study clubs and societies to specifically meet the educational needs of practicing dentists.



www.gidedental.com

gIDE by the numbers:

108+
Days of
classroom
training

40,000+
Newsletter
subscribers

25,000+
Monthly
visitors

10,000+
Registered
dentists

2,200+
Master clinician
program
graduates

100+
Countries of
members

For more information visit **HuFriedyGroup.com**