

POWER



AIR-FLOW® HANDY 3.0 DEVICES


THE BETTER WAY FOR
BIOFILM MANAGEMENT

Introducing **AIR-FLOW[®]** **THERAPY**

AIR-FLOW Therapy is completely different than other air polishing systems. This advanced technology offers benefits to both the clinician and patient that are difficult to replicate with other procedures. AIR-FLOW Therapy, along with glycine PERIO powder, is safe and effective for tissue, enamel, root surfaces, porcelain composite restorations, orthodontic appliances and implants & implant abutments.



WHAT IT DOES:




- Polishes teeth
- Removes stain
- Disrupts subgingival biofilm

PATIENT BENEFITS:

- Increased comfort
- Healthy gingiva
- Grit-free polishing
- Faster procedure

CLINICIAN BENEFITS:

- Effective biofilm disruption
- Efficient treatments
- Happier patients



BRING AIR-FLOW TO YOUR PRACTICE WITH A HANDY DEVICE! A great introduction to AIR-FLOW Therapy, the handy 3.0 is a lightweight handheld device that connects to your dental unit. Designed with a team of ergonomics experts, the slim handpiece is easy for the clinician to grip and rotate, while the position of the powder chamber ensures an unrestricted view of the treatment area.

CLINICAL CAPABILITIES



STANDARD HANDPIECE

- Designed for Classic and Classic Comfort powders
- Ideal for supragingival stain removal applications
- Included with standard AIR-FLOW handy 3.0 device



PLUS HANDPIECE

- Designed for fine-grain PERIO powder
- Ideal for supragingival biofilm removal, including on ortho and restorative materials, as well as subgingival biofilm removal up to 4mm
- Safe for perio and implant maintenance patients
- Included with AIR-FLOW handy 3.0 PREMIUM device and compatible with AIR-FLOW handy 3.0 PERIO device (sold separately)



PERIO HANDPIECE

- Designed for fine-grain PERIO Powder
- Flexible plastic nozzle is inserted directly into pockets to flush out biofilm
- Safe for perio and implant maintenance patients
- Included with AIR-FLOW handy 3.0 PREMIUM and AIR-FLOW handy 3.0 PERIO devices


DEVICE OPTIONS



AIR-FLOW HANDY 3.0 PREMIUM

- For complete biofilm removal
- Includes both PLUS and PERIO handpieces for supragingival and subgingival applications
- Easily switch between handpieces during treatment
- Designed for optimal delivery of PERIO Powder

NEW



AIR-FLOW HANDY 3.0 PERIO

- For direct subgingival biofilm removal
- Includes PERIO handpiece for subgingival applications up to 5mm deep*
- Designed for PERIO Powder

*AIR-FLOW handy 3.0 Perio is approved for use in pockets up to 10mm in Canada



AIR-FLOW HANDY 3.0

- For heavy stain removal
- Includes standard handpiece for mainly supragingival applications
- Designed for Classic and Classic Comfort Powders, but can be used with PERIO Powder



PERIO POWDER

Perio powder's extra-fine, low-density grains safely and completely remove harmful biofilm and bacteria both supragingivally and subgingivally. Made from glycine, PERIO powder is gentle on enamel and restorative surfaces, making the patient experience more comfortable.

COMFORT POWDER

The Classic Comfort powder offers the well-known polishing efficiency of sodium bicarbonate with increased patient comfort due to a smaller than average grain size.

While the Classic Comfort powder is primarily used for supragingival removal of biofilm and stains, it can also be used for surface preparation before dental treatments.

CLASSIC POWDER

The AIR-FLOW Classic powder is effective for dental cleanings and plaque removal. Made from sodium bicarbonate, the most commonly used material for air polishing, AIR-FLOW Classic powders provide an alkaline effect to protect against tooth decay, and an osmotic effect to support the treatment of gingivitis.

Available flavors: Lemon and Mint

ORDER YOUR HANDY DEVICE!

PRODUCT	PART CODES*
AIR-FLOW handy 3.0 Midwest*	FT-220#HF/A/001
AIR-FLOW handy 3.0 Kavo	FT-220#HF/A/002
AIR-FLOW handy 3.0 PERIO Midwest*	FT-221#HF/A/001
AIR-FLOW handy 3.0 PERIO Kavo	FT-221#HF/A/002
AIR-FLOW handy 3.0 Premium Midwest*	FT-225#HF/A/001
AIR-FLOW handy 3.0 Premium Kavo	FT-225#HF/A/002
handy 3.0 Handpiece	EL-540#HF/A
handy 3.0 PERIO Handpiece	EL-542#HF/A
handy 3.0 PLUS Handpiece	EL-545#HF/A

NOZZLE REFILLS	PART CODES
PERIO-FLOW nozzle refill box (40ct)	DT-113HF
PERIO-FLOW nozzle refill pack (4 boxes of 40)	DT-113HF/HF


POWDERS IN PACKS OF 4	PART CODES
AIR-FLOW Classic Powder, Lemon	DV-111/A/LEM
AIR-FLOW Classic Powder, Mint	DV-111/A/MIN
AIR-FLOW Classic Comfort Powder	DV-109/A
AIR-FLOW PERIO Powder	DV-110/A/HS

*Fits standard, 4-hole connection

*W&H, NSK, Sirona and Bien Air Connections also available

For more information, visit

Hu-Friedy.com/BIOFILM



Hu-Friedy Mfg. Co., LLC: 3232 N. Rockwell St. | Chicago, IL 60618 | USA
 Telephone: 1-800-Hu-Friedy (1-773-975-6100) | Website: Hu-Friedy.com

©2018 Hu-Friedy Mfg. Co., LLC. All rights reserved. [557] HF-3030/0118

NOTE: Kavo, Midwest, W&H, NSK, Bien Air and Sirona are trademarks or registered trademarks of their respective owners.

